

Product Datasheet

Rat Osm protein (orb605273)

Description	Recombinant Rat Oncostatin-M(Osm)
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 10xHis-tagged
Note	For research use only
Application notes	Full Length of Mature Protein
Protein Sequence	KRGCSSSSPKLLSQLKSQANITGNTASLLEPYILHQNLNLTLLRAACTEHPVAFPSDMLRQL SKPDFLSTVHATLGRVWHQLGAFRQQFPKIQDFPELERARQNIQGIRNNVYCMARLLHPPL EIPEPTQADSGTSRPTTTAPGIFQIKIDSCRFLWGYHRFMGSVGRVFEEWGDGSRRSRR
Purity	Greater than 85% as determined by SDS-PAGE.
Source	Mammalian cell
MW	24.2 kDa
Uniprot ID	Q65Z15
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)