
Product Datasheet

PSAP Antibody / Prostate Specific Acid Phosphatase (orb606367)

Description

Prostate specific acid phosphatase (PSAP) is an

orthophospho...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

IHC-P

Immunogen

Human recombinant full length protein was used as the immunogen for this recombinant PSAP antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

The recombinant PSAP antibody (with azide) can be stored at 2-8°C. The azide-free format should be aliquoted and stored at -20°C or colder.

Note

For research use only

Application notes

Titering of the recombinant PSAP antibody may be required for optimal performance.1. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Recombinant

Purity

Protein G affinity chromatography

Uniprot ID

P15309

Hazard Information

This recombinant PSAP antibody is available for research use only.

Entrez


55

Dilution Range

Immunohistochemistry (FFPE): 0.5-1ug/ml for 30 min at RT


Expiration Date

12 months from date of receipt.



IHC testing of human prostate carcinoma

...



SDS-PAGE analysis of purified, BSA-free ...