
Product Datasheet

KRT19 Antibody / Cytokeratin 19 (orb606439)

Description

Recognizes a protein of 40kDa, identified as cytokeratin-


Species/Host

Mouse

Reactivity


Human

Conjugation

Unconjugated

Tested Applications

FACS, IHC-P, WB



IHC testing of FFPE human colon tissue w...

Immunogen


Recombinant human KRT19 protein was used as the immunogen for the recombinant Cytokeratin 19 antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the recombinant Cytokeratin 19 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).



Flow cytometry testing of PFA-fixed huma...

Note

For research use only

Application notes

Optimal dilution of the recombinant Cytokeratin 19 antibody should be determined by the researcher.1. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Recombinant

Purity

Protein G affinity chromatography

Uniprot ID

P08727

Hazard Information

This recombinant Cytokeratin 19 antibody is available for research use only.

Dilution Range

Flow cytometry: 1-2ug/million cells,Western blot: 1-2ug/ml,Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT



Western blot testing of human HepG2 cell...

Expiration Date

12 months from date of receipt.