

Product Datasheet

CD43 Antibody (orb606475)

Description

It recognizes a cell surface glycoprotein of


Species/Host

Mouse

Reactivity


Human

Conjugation

Unconjugated

Tested

FACS, IF, IHC-P, WB

Applications


IHC testing of FFPE human tonsil with re...

Immunogen

Recombinant human protein was used as the immunogen for the recombinant CD43 antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage


Store the recombinant CD43 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the recombinant CD43 antibody should be determined by the researcher.1. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.



IHC testing of FFPE human lymphoma with ...

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Recombinant

Purity

Protein G affinity chromatography

Uniprot ID

P16150

Hazard
Information


This recombinant CD43 antibody is available for research use only.

Dilution Range

Western blot: 1-2ug/ml, Flow cytometry: 1-2ug/million cells, Immunofluorescence: 1-2ug/ml, Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

Expiration Date

12 months from date of receipt.



Immunofluorescence staining of human K56...