

---

## Product Datasheet

### NKX2.1 Antibody / TTF1 (orb606478)

**Description**

NKX2.1/TTF-1 is a member of the NKx2 family of


**Species/Host**

Rabbit

**Reactivity**

Human, Mouse, Rat

**Conjugation**

Unconjugated

**Tested Applications**

IHC-P

**Immunogen**

Recombinant protein was used as the immunogen for the recombinant NKX2.1 antibody.

**Preservatives**

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

**Storage**


Store the recombinant NKX2.1 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

**Note**

For research use only

**Application notes**

Optimal dilution of the recombinant NKX2.1 antibody should be determined by the researcher.1. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.



IHC:  
Formalin-fixed, paraffin-embedded h...

**Formula**

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

**Isotype**

Rabbit IgG, kappa

**Clonality**

Recombinant

**Purity**

Protein A affinity chromatography

**Uniprot ID**

**P43699**

**Hazard Information**


This recombinant NKX2.1 antibody is available for research use only.

**Dilution Range**


Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

**Expiration Date**

12 months from date of receipt.



IHC:  
Formalin-fixed, paraffin-embedded h...



SDS-PAGE analysis of purified, BSA-free ...