
Product Datasheet

XAB1 Antibody (orb606573)

Description

GPN1 (GPN-loop GTPase 1), also called XAB1 (XPA-binding


Species/Host

Mouse

Reactivity


Human

Conjugation

Unconjugated

Tested Applications

ELISA, WB



IHC testing of FFPE human pancreas with

...

Immunogen

Recombinant full length human protein was used as the immunogen for this XAB1 antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the XAB1 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).


Note

For research use only

Application notes

The concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the XAB1 antibody to be titered up or down for optimal performance.

Human Protein Microarray Specificity Validation



Analysis of HuProt (TM) microarray conta...

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Purity


Protein G affinity chromatography

Uniprot ID

[Q9HCN4](#)

Hazard Information

This XAB1 antibody is available for research use only.



SDS-PAGE analysis of purified, BSA-free ...

Dilution Range

ELISA (order BSA/sodium azide-free format for coating), Western blot: 1-2ug/ml

Expiration Date

12 months from date of receipt.