



## Product Datasheet

### Aurora B Antibody / AURKB (orb606677)


<b>Description</b>	Recognizes a protein of 39kDa, which is identified as Aurora B. The serine/threonine protein kinase ...
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC-P
<b>Immunogen</b>	Amino acids 89-251 were used as the immunogen for this Aurora B antibody.
<b>Preservatives</b>	0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide
<b>Storage</b>	Store the Aurora B antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).
<b>Note</b>	For research use only
<b>Application notes</b>	Optimal dilution of the Aurora B antibody should be determined by the researcher.
<b>Formula</b>	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide
<b>Isotype</b>	Mouse IgG1, kappa
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Protein G affinity
<b>Uniprot ID</b>	<a href="#">Q96GD4</a>
<b>Hazard Information</b>	This Aurora B antibody is available for research use only.
<b>Dilution Range</b>	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
<b>Expiration Date</b>	12 months from date of receipt.



IHC testing of FFPE human prostate tissue...



SDS-PAGE analysis of purified, BSA-free ...



Analysis of HuProt (TM) microarray conta...