
Product Datasheet

TRP1 / Tyrosinase-Related Protein-1 Antibody (orb606789)

Description

Tyrosinase-Related Protein-1 / TYRP1 / TRP1 is

m...

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

IHC-P

Immunogen

Human recombinant protein used as the immunogen for this recombinant TYRP1 antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the recombinant TYRP1 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

The optimal dilution of the recombinant TYRP1 antibody for each application should be determined by the researcher.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Rabbit IgG, kappa

Clonality

Recombinant

Purity

Protein A affinity chromatography

Uniprot ID
P17643
Hazard Information


This recombinant TYRP1 antibody is available for research use only.

Dilution Range


Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

Expiration Date

12 months from date of receipt.



IHC testing of FFPE human melanoma with ...



SDS-PAGE analysis of purified, BSA-free ...