

Product Datasheet

**TTR Antibody / Transthyretin / Prealbumin
(orb606853)**

Description

Prealbumin, also designated transthyretin, is a

Human Protein Microarray Specificity Validation



	tr...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant human protein was used as the immunogen for this Prealbumin antibody.
Preservatives	0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide
Storage	Store the Prealbumin antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).
Note	For research use only
Application notes	Optimal dilution of the Prealbumin antibody should be determined by the researcher.
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide
Isotype	Mouse IgG2a, kappa
Clonality	Monoclonal
Purity	Protein G affinity chromatography
Uniprot ID	P02766
Hazard Information	This Prealbumin antibody is available for research use only.
Dilution Range	ELISA (order BSA/sodium azide-free format for coating),Western blot: 1-2ug/ml
Expiration Date	12 months from date of receipt.

Analysis of HuProt (TM) microarray conta...



Western blot testing of human liver lysa...



Western blot testing of human kidney lys...