

## Product Datasheet

### SULT1E1 Antibody / Estrogen Sulfotransferase (orb606978)

**Description**

The soluble sulfotransferases contribute to the


**Species/Host**

Mouse

**Reactivity**

Human

**Conjugation**

Unconjugated

**Tested**

IHC-P, WB

**Applications**

Western  
blot testing  
of human  
liver lysa...

**Immunogen**

Recombinant human protein was used as the immunogen for the SULT1E1 antibody.

**Preservatives**

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

**Storage**

Store the SULT1E1 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

**Note**

For research use only

**Application notes**

Optimal dilution of the SULT1E1 antibody should be determined by the researcher.1. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.

IHC staining  
of FFPE  
human  
small  
intesti...

**Formula**

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

**Isotype**

Mouse IgG2b, kappa

**Clonality**

Monoclonal

**Purity**

Protein G affinity chromatography

**Uniprot ID**

**P49888**

**Hazard**
**Information**


This SULT1E1 antibody is available for research use only.

**Dilution Range**

Immunohistochemistry (FFPE): 0.5-1ug/ml for 30 min at RT, Western blot: 1-2ug/ml

**Expiration Date**

12 months from date of receipt.



IHC staining  
of FFPE  
human  
lymph node  
wi...