
Product Datasheet

Ku (phospho-70/80) antibody (CF488) (orb607366)

Description

Mouse monoclonal antibody against Ku

CF488

Species/Host

Mouse

Reactivity

Human, Primate, Rat

Conjugation

Unconjugated

Tested Applications

FC, IF

Immunogen

Nuclear extract of human HL-60 cells

Target

XRCC6 (p70) & XRCC5 (p80)

Preservatives

Antibody Purified from Bioreactor Concentrate by Protein A/G and conjugated to various reporter molecules. Prepared in 10mM PBS with 0.05% rAlbumin and 0.05% azide. Contact us if you require this Ab in a different format.

Concentration

Purified Ab conjugated to CF488

Storage

Antibody with azide - store at 4 to 8°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Isotype

Mouse / IgG1, kappa

Clonality

Monoclonal

MW

70kDa & 80kDa

Uniprot ID
[P13010](#), [P12956](#)
Entrez
[7520](#), [2547](#)
Dilution Range

Flow Cytometry (5ul per test per one million cells or 5ul per 100ul of whole blood);Immunofluorescence (1:50-1:100);Optimal dilution for a specific application should be determined.

Expiration Date

12 months from date of receipt.