

Product Datasheet

S100A9 + Calprotectin antibody (orb607694)

Description

Rabbit monoclonal antibody to S100A9 +

Species/Host	Rabbit
Reactivity	Bovine, Canine, Equine, Feline, Goat, Guinea pig, Human, Monkey, Mouse, Porcine, Primate, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Affinity Purified monocyte membrane preparation
Target	S100A8 & S100A9
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Isotype	Rabbit / IgG
Clonality	Monoclonal
MW	12-14kDa (doublet)
Uniprot ID	P06702, P05109
Entrez	6280, 6279
Dilution Range	Immunohistochemistry (Formalin-fixed) (1-2ug/ml for 30 minutes at RT)(Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes); Optimal dilution for a specific application should be determined.
Expiration Date	12 months from date of receipt.