
Product Datasheet

CLEC9A antibody (orb607928)

Description

Rat monoclonal antibody to CLEC9A

Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	RBL-2H3 cells expressing Mouse CLEC9A fused to an HA epitope.
Target	CLEC9A (Mouse)
Preservatives	200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Isotype	Rat / IgG1
Clonality	Monoclonal
MW	27kDa
Uniprot ID	Q8BRU4
Entrez	232414
Dilution Range	Flow Cytometry (0.5-1ug/million cells); ,Western Blot (1-2ug/ml);,Immunofluorescence (0.5-1ug/ml); ,Optimal dilution for a specific application should be determined.
Expiration Date	12 months from date of receipt.