
Product Datasheet

Pseudomonas aeruginosa serotype 5C antibody (orb607942)

Description

Mouse monoclonal antibody to

Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Pseudomonas aeruginosa serotype 5C
Target	Not Applicable
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Isotype	Mouse / IgM, kappa
Clonality	Monoclonal
MW	Not Known
Uniprot ID	Not Applicable
Entrez	Not Applicable
Dilution Range	ELISA (For coating use Ab at 1-5ug/ml, order Ab without BSA);Immunofluorescence (1-2ug/ml); Optimal dilution for a specific application should be determined.
Expiration Date	12 months from date of receipt.