

---

## Product Datasheet

### HIV1-p24 antibody (CF488) (orb607953)

**Description**

Mouse monoclonal antibody against Human

p24) conjugated to C...

**Species/Host**

Mouse

**Conjugation**

Unconjugated

**Tested Applications**

ICC, IF

**Immunogen**

Recombinant HIV-1 Gag p24 protein

**Preservatives**

Antibody Purified from Bioreactor Concentrate by Protein A/G and conjugated to various reporter molecules. Prepared in 10mM PBS with 0.05% rAlbumin and 0.05% azide. Contact us if you require this Ab in a different format.

**Concentration**

Purified Ab conjugated to CF488

**Storage**

Antibody with azide - store at 4 to 8°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

**Note**

For research use only

**Isotype**

Mouse / IgG1, kappa

**Clonality**

Monoclonal

**MW**

24kDa (mature); 55kDa & 41kDa (precursors)

**Uniprot ID**

**P04585**

**Entrez**

**155348**

**Dilution Range**

Immunocytochemistry (Acetone or paraformaldehyde fixed) (1:50-1:100);Immunofluorescence (1:50-1:100);Optimal dilution for a specific application should be determined.

**Expiration Date**

12 months from date of receipt.