

---

## Product Datasheet

### HPV-16 antibody (CF488) (orb607958)

**Description**

Mouse monoclonal antibody against HPV-16

	CF488
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FC, IF
<b>Immunogen</b>	Human papilloma virus type 16, major capsid protein L1
<b>Target</b>	Not Applicable
<b>Preservatives</b>	Antibody Purified from Bioreactor Concentrate by Protein A/G and conjugated to various reporter molecules. Prepared in 10mM PBS with 0.05% rAlbumin and 0.05% azide. Contact us if you require this Ab in a different format.
<b>Concentration</b>	Purified Ab conjugated to CF488
<b>Storage</b>	Antibody with azide - store at 4 to 8°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
<b>Note</b>	For research use only
<b>Isotype</b>	Mouse / IgG2a, kappa
<b>Clonality</b>	Monoclonal
<b>MW</b>	57kDa
<b>Uniprot ID</b>	Not Applicable
<b>Entrez</b>	Not Applicable
<b>Dilution Range</b>	Flow Cytometry (5ul per test per one million cells or 5ul per 100ul of whole blood);Immunofluorescence (1:50-1:100);Optimal dilution for a specific application should be determined.
<b>Expiration Date</b>	12 months from date of receipt.