
Product Datasheet

LMP-1 antibody (orb607959)

Description

Mouse monoclonal antibody to Epstein-Barr

Species/Host

Mouse

Conjugation

Unconjugated

Tested

ELISA, WB

Applications**Immunogen**

Recombinant fusion protein containing the sequence of bacterial beta-galactosidase and the carboxyl half of EBV-encoded LMP

Target

Not Applicable

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Isotype

Mouse / IgG1, kappa

Clonality

Monoclonal

MW

60kDa

Uniprot ID

Not Applicable

Entrez

Not Applicable

Dilution Range

ELISA (For coating, purchase antibody without BSA); Western Blot (1-2ug/ml); Optimal dilution for a specific application should be determined.

Expiration Date

12 months from date of receipt.