
Product Datasheet

HPV-16 E6 + HPV-18 E6 antibody (orb607997)

Description

Mouse monoclonal antibody to HPV-16

16/18)

Species/Host

Mouse

Conjugation

Unconjugated

Tested Applications

ELISA

Immunogen

 HPV18 E6- β -galactosidase fusion protein

Target

Not Applicable

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Isotype

Mouse / IgG1, kappa

Clonality

Monoclonal

MW

16/17kDa

Uniprot ID

Not Applicable

Entrez

Not Applicable

Dilution Range

ELISA (For coating, order antibody without BSA); Optimal dilution for a specific application should be determined.

Expiration Date

12 months from date of receipt.