

Product Datasheet

ATM (phospho-Ser1981) antibody (orb608078)

Description	ATM (phospho-Ser1981) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human ATM around the phosphorylation site of Ser1981 (EG(p-S)Q)
Target	phospho-ATM (Ser1981)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	370 kDa
Uniprot ID	Q62388
Dilution Range	WB=1:300-1000, IHC-P=1:50-100, ICC=1:50-100, IHC-F=1:50-100, IF=1:50-100
Expiration Date	12 months from date of receipt.