

Product Datasheet

Liraglutide (Victoza) ELISA Kit (orb608237)

Description	Liraglutide (Victoza) ELISA Kit
Reactivity	Human
Range	40 ng/ml - 2560 ng/ml
Tested Applications	ELISA
Form/Appearance	This kit comes with a pre-coated, standard 96 well plate, along with relevant detection and reaction reagents (TMB Substrate, Stop Solution).
Note	For research use only
Application notes	Sandwich Assay
Sample Types	Human serum and plasma
Sensitivity	less than 40ng/ml
Expiration Date	Please enquire.