

## Product Datasheet

### Lixisenatide ELISA Kit (orb608238)

<b>Description</b>	Lixisenatide ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.1-1000 ng/ml
<b>Tested Applications</b>	ELISA
<b>Form/Appearance</b>	This kit comes with a pre-coated, standard 96 well plate, along with relevant detection and reaction reagents (TMB Substrate, Stop Solution).
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich Assay
<b>Sample Types</b>	Human serum, plasma and tissue extracts
<b>Sensitivity</b>	~180 ng/ml.
<b>Expiration Date</b>	Please enquire.