



Product Datasheet

TLR4 antibody (orb609129)


Description	Goat polyclonal antibody to TLR4
Species/Host	Goat
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC
Target	TLR4 / CD284
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Protein Sequence	YSSQDEDWVRNE
Clonality	Polyclonal
MW	95.7
Entrez	7099
Dilution Range	P-ELISA: 1:16000, IF/ICC: 10 µg/ml
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of NIH3T3 ce...



orb609129 Immunofluorescence analysis of...



orb609129 Flow cytometric analysis of pa...