

## Product Datasheet

### Lisinopril Dihydrate (orb609430)

<b>Description</b>	Lisinopril Dihydrate
<b>Storage</b>	Storage Temp: Ambient
<b>Note</b>	For research use only
<b>Formula</b>	$C_{21}H_{31}N_3O_5 \cdot 2H_2O$
<b>Purity</b>	≥98%
<b>MW</b>	441.52
<b>CAS Number</b>	83915-83-7
<b>SMILES</b>	<chem>C1CC(N(C1)C(=O)C(CCCCN)NC(CCC2=CC=CC=C2)C(=O)O)C(=O)O.O.O</chem>
<b>Solubility (25°C)</b>	Soluble in water (100 mg/mL) or DMSO (6.5mg/mL with heating and sonication), methanol (14 mg/mL). Insoluble in ethanol.
<b>Hazard Information</b>	Chemical Hazard Code: None. UN Number: Not dangerous goods
<b>Expiration Date</b>	12 months from date of receipt.