

Product Datasheet

Tetrahydroberberine (orb609704)

Description	Tetrahydroberberine
Storage	Storage Temp: 4°C
Note	For research use only
Formula	$C_{20}H_{21}NO_4$
Purity	≥98%
MW	339.39
CAS Number	522-97-4
SMILES	<chem>COc1ccc2CC3N(CCc4cc5OCOc5cc34)Cc2c1OC</chem>
Solubility (25°C)	10 mM in DMSO insoluble in water
Hazard Information	UN Number: Not dangerous goods
Expiration Date	12 months from date of receipt.