

Product Datasheet

TG100-115 (orb609725)

Description	TG100-115
Storage	Storage Temp: -20°C
Note	For research use only
Formula	$C_{18}H_{14}N_6O_2$
Purity	≥98%
MW	346.34
CAS Number	677297-51-7
SMILES	<chem>c1cc(cc(c1)O)c2c(nc3c(n2)c(nc(n3)N)N)c4cccc(c4)O</chem>
Solubility (25°C)	DMSO 9 mg/mL (25.98 mM) Water Insoluble Ethanol Insoluble
Hazard Information	UN Number: Not dangerous goods
Expiration Date	12 months from date of receipt.