

Product Datasheet

SGX-523 (orb609816)

Description	SGX-523
Storage	Storage Temp: -20°C
Note	For research use only
Formula	$C_{18}H_{13}N_7S$
Purity	≥98%
MW	359.41
CAS Number	1022150-57-7
SMILES	<chem>Cn1cc(cn1)c2ccc3nnc(n3n2)Sc4ccc5c(c4)cccn5</chem>
Solubility (25°C)	DMSO 3 mg/mL warmed (8.34 mM) Water Insoluble Ethanol Insoluble
Hazard Information	UN Number: Not dangerous goods
Expiration Date	12 months from date of receipt.