

Product Datasheet

AdipoRon (orb610097)

Description	AdipoRon
Storage	Storage Temp: -20°C
Note	For research use only
Purity	≥98%
MW	428.53
CAS Number	924416-43-3
SMILES	<chem>c1ccc(cc1)CN2CCC(CC2)NC(=O)COc3ccc(cc3)C(=O)c4ccccc4</chem>
Solubility (25°C)	Soluble in DMSO, 100% ethanol or methanol. Insoluble in water
Hazard Information	UN Number: Not dangerous goods
Expiration Date	12 months from date of receipt.