



Product Datasheet

mtTFA TFAM Antibody (monoclonal, 4D9) (orb610917)


Description	mtTFA TFAM Antibody (monoclonal, 4D9)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human mtTFA (214-241aa EMKSWEEQMIEVGRKDLLRRTIKKQRKY), different from the related mouse and rat sequences by five amino acids.
Form/Appearance	Lyophilized
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Add 0.2ml of distilled water will yield a concentration of 500µg/ml.
Isotype	Mouse IgG2b
Clonality	Monoclonal
MW	40 kDa
Uniprot ID	Q00059
Expiration Date	12 months from date of receipt.



Flow Cytometry analysis of CACO-2 cells ...



IHC analysis of TFAM using anti-TFAM ant...



IHC analysis of TFAM using anti-TFAM ant...