

Product Datasheet

HIV1 p24 antibody (orb610960)

Description	Mouse monoclonal antibody to HIV1 p24
Reactivity	Virus
Conjugation	Unconjugated
Immunogen	r.HIV-1p24 Core
Preservatives	Lyophilized from PBS.
Form/Appearance	Sterile Filtered White lyophilized (freeze-dried) powder.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Monoclonal Antibodies, Anti-Viral
Isotype	mouse IgG1
Clonality	Monoclonal
Solubility (25°C)	Reconstitute with 1ml 1x PBS. Mix gently, wash the sides of the vial and wait 30-60 seconds before use.
Expiration Date	12 months from date of receipt.