

## Product Datasheet

### Human IL34 Protein (orb610989)

<b>Description</b>	Recombinant Human IL34 protein
<b>Preservatives</b>	IL34 protein solution (0.25mg/ml) contains 20mM Tris-Hcl buffer (pH8.5) 1mM DTT, 0.15M NaCl, 0.1mM PMSF, 1mM EDTA and 30% glycerol.
<b>Form/Appearance</b>	Sterile Filtered colorless solution.
<b>Storage</b>	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
<b>Note</b>	For research use only
<b>Application notes</b>	Cytokines And Growth Factors
<b>Protein Sequence</b>	ADPNEPLEMW PLTQNEECTV TGFLRDKLQY RSRLQYMKHY FPINYKISVP YEGVFRIANV TRLQRAQVSE RELRYLWVLV SLSATESVQD VLEGHPSWK YLQEVQTL L NVQQGLTDVE VSPKVESVLS LLNAPGPNLK LVRPKALLDN CFRVMELLYC SCCKQSSVLN WQDCEVPSPO SCSPEPSLQY AATQLYPPPP WSPSSPPHST GSVRPVRAQG EGLLPHHHHH H
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE.
<b>Source</b>	Sf9, Insect cells
<b>Expiration Date</b>	12 months from date of receipt.