

Product Datasheet

Rat IL1 alpha Protein (orb611017)

Description	Recombinant Rat IL1 alpha protein
Preservatives	IL1A protein solution (1mg/ml) contains 10% glycerol & Phosphate Buffered Saline (pH 7.4).
Form/Appearance	Sterile Filtered colorless solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Cytokines And Growth Factors
Protein Sequence	MGSSHHHHHH SSGLVPRGSH MGSSAPHSFQ NNLRYKLIRI VKQEFIMNDS LNQNIYVDMMDRIHLKAASLN DLQLEVKFDM YAYSSGGDDSKYPVTLKVSNTQLFVSAQGE DKPVLLKEIPETPKLITGSE TDLIFFWEKI NSKNYFTSAA FPELLIATKE QSQVHLARGL PSMIDFQIS
Purity	Greater than 95% as determined by SDS-PAGE.
Source	E.coli
Expiration Date	12 months from date of receipt.