

## Product Datasheet

### OSI-027 (orb611257)

|                        |   |
|------------------------|---|
| <b>Description</b>     | OSI-027 is a selective and potent dual inhibitor of mTORC1 and mTORC2 with IC50 of 22 nM and 65 nM, ... |
| <b>Note</b>            | For research use only   |
| <b>Formula</b>         | $C_{21}H_{22}N_6O_3$  |
| <b>MW</b>              | 406.44  |
| <b>CAS Number</b>      | [936890-98-1]   |
| <b>SMILES</b>          | <chem>O=C([C@H]1CC[C@H](C2=NC(C(N3)=CC4=C3C(OC)=CC=C4)=C5C(N)=NC=NN52)CC1)O</chem>                      |
| <b>Expiration Date</b> | 12 months from date of receipt.   |