

Product Datasheet

YKL-05-099 (orb611258)

Description	YKL-05-099 is a salt-inducible kinase (SIK) probe; inhibits SIK2 with an IC50 of 40 nM.
Note	For research use only
Formula	$C_{32}H_{34}ClN_7O_3$
MW	600.11
CAS Number	[1936529-65-5]
SMILES	<chem>O=C1N(C2=NC=C(OC)C=C2)C3=NC(NC4=CC=C(C5CCN(C)CC5)C=C4O)C=NC=C3CN1C6=C(C)C=CC=C6Cl</chem>
Expiration Date	12 months from date of receipt.