

## Product Datasheet

### SEC inhibitor KL-1 (orb611271)

<b>Description</b>	SEC inhibitor KL-1 is a peptidomimetic inhibitor of super elongation complex (SEC) and transcription...
<b>Note</b>	For research use only
<b>Formula</b>	$C_{18}H_{16}ClNO_4$
<b>MW</b>	345.08
<b>CAS Number</b>	[900308-84-1]
<b>SMILES</b>	<chem>C1=C(OC)C=CC=C1C(CC(C(=O)NC1=CC(Cl)=CC=C1C)=O)=O</chem>
<b>Expiration Date</b>	12 months from date of receipt.