

Product Datasheet

CC-671 (orb611300)

Description	CC-671 is a potent, selective, dual TTK (Mps1)/CLK2 inhibitor with IC50 of 5/3 nM, respectively.
Note	For research use only
Formula	$C_{28}H_{28}N_6O_4$
MW	512.57
CAS Number	[1618658-88-0]
SMILES	<chem>CC1=NC2=C(O1)C=C(C=C2)C3=CNC4=C3C(=NC(=N4)NC5=C(C=C(C=C5)C(=O)NC)OC)OC6CCCC6</chem>
Expiration Date	12 months from date of receipt.