

Product Datasheet

Bobcat339 (orb611313)

Description	Bobcat339 (Bobcat 339) is a novel cytosine-based TET enzyme inhibitor with IC50 of 33 uM (TET1) and ...
Note	For research use only
Formula	$C_{16}H_{12}ClN_3O$
MW	297.7389822
CAS Number	[2280037-51-4]
SMILES	<chem>C1C(N([H])[H])=NC(N(C=1[H])C1=C([H])C([H])=C([H])C(C2C([H])=C([H])C([H])=C([H])C=2[H])=C1[H])=O</chem>
Expiration Date	12 months from date of receipt.