

Product Datasheet

E3 ligase Ligand 1A (orb611322)

Description	E3 ligase Ligand 1A is a ligand of E3 ligase, used in PROTAC technology; E3 ligase Ligand 1A can be ...
Note	For research use only
Formula	$C_{27}H_{28}N_2O_5S_2$
MW	444.59
CAS Number	[1948273-02-6]
SMILES	<chem>N[C@@H](C(C)(C)C)C(N1[C@H](C(N[C@@H](C)C2=CC=C(C3=C(C)N=C(S3)C=C2)=O)C[C@@H](O)C1)=O</chem>
Expiration Date	12 months from date of receipt.