

Product Datasheet

DG051 (orb611348)

Description	A potent, selective LTA4H (Leukotriene A4 hydrolase) inhibitor (IC50/Kd=47/25 nM) of leukotriene B4 ...
Note	For research use only
Formula	$C_{21}H_{25}Cl_2NO_4$
MW	426.3335
CAS Number	[929915-58-2]
SMILES	<chem>C1C[C@H](N(C1)CCCC(=O)O)COC2=CC=C(C=C2)OC3=CC=C(C=C3)Cl.</chem> Cl
Expiration Date	12 months from date of receipt.