

## Product Datasheet

### MP-A08 (orb611453)

<b>Description</b>	MP-A08 is a first-in-class, highly selective, ATP competitive sphingosine kinase (SphK) inhibitor (K...
<b>Note</b>	For research use only
<b>Formula</b>	$C_{27}H_{25}N_3O_4S_2$
<b>MW</b>	519.6351
<b>CAS Number</b>	[219832-49-2]
<b>SMILES</b>	<chem>O=S(=O)(NC1C(C=NC2C(NS(=O)(C3C=CC(C)=CC=3)=O)=CC=CC=2)=CC=CC=1)C1C=CC(C)=CC=1</chem>
<b>Expiration Date</b>	12 months from date of receipt.