

Product Datasheet

LUF6000 (orb611495)

Description	LUF6000 is a potent, selective, positive allosteric modulator (enhancer) of human A3 adenosine recep...
Note	For research use only
Formula	$C_{22}H_{20}Cl_2N_4$
MW	411.327
CAS Number	[890087-21-5]
SMILES	<chem>C1C(Cl)=CC=C(NC2C3=C(N=C(C4CCCC4)N3)C3=C(C=CC=C3)N=2)C=1</chem>
Expiration Date	12 months from date of receipt.