

Product Datasheet

MK-8033 (orb611565)

| | |
|------------------------|---|
| Description | MK-8033 is a potent, specific, dual c-Met/Ron kinases inhibitor with IC50 of 1/7 nM, also potently i... |
| Note | For research use only |
| Formula | $C_{25}H_{21}N_5O_3S$ |
| MW | 471.5309436 |
| CAS Number | [1001917-37-8] |
| SMILES | <chem>CN1C=C(C2=CN=C3C=CC4C=CC(CS(NCC5C=CC=CN=5))(=O)=O)=CC=4C(=O)C3=C2)C=N1</chem> |
| Expiration Date | 12 months from date of receipt. |