

Product Datasheet

Zaprinast (orb611672)

Description	Zaprinast (M&B 22948) is a potent, cGMP-specific phosphodiesterase inhibitor with IC50 of 0.5-0.76 a...
Note	For research use only
Formula	$C_{13}H_{13}N_5O_2$
MW	271.2746217
CAS Number	[37762-06-4]
SMILES	<chem>CCCOC1=CC=CC=C1C1=NC(=O)C2=NNNC2=N1</chem>
Expiration Date	12 months from date of receipt.