

Product Datasheet

Bay K-8644 (orb611678)

Description	Bay K-8644 is a potent L-type calcium channels (LTCCs) activator with EC50 of 17.3 nM, has positive ...
Note	For research use only
Formula	$C_{16}H_{15}N_2O_4F_3$
MW	356.2965
CAS Number	[71145-03-4]
SMILES	<chem>CC1=C(C(C2=CC=CC=C2C(F)(F)F)C(=C(C)N1)[N+](=O)[O-])C(=O)OC</chem>
Expiration Date	12 months from date of receipt.