

Product Datasheet

CGP-42112 (orb611698)

Description	CGP-42112 (CGP42112A) is a highly potent, selective Angiotensin-II subtype 2 receptor(AT2 receptor) ...
Note	For research use only
Formula	$C_{52}H_{69}N_{13}O_{11}$
MW	1052.185
CAS Number	[127060-75-7]
SMILES	<chem>N/C(=N/CCC[C@@H](C(NCCCC[C@@H](C(N[C@H](C(N1CCC[C@H]1C(N[C@@H]([C@H](CC)C(=O)O)=O)=O)CC1=CN=CN1)=O)NC([C@@H](NC(C1C=CC=NC=1)=O)CC1C=CC(O)=CC=1)=O)=O)NC(OCC1C=CC=CC=1)=O)/N</chem>
Expiration Date	12 months from date of receipt.