

## Product Datasheet

### JAK3-IN-1 (orb611792)

<b>Description</b>	JAK3-IN-1 is a potent and selective JAK3 inhibitor with IC50 of 4.8 nM.
<b>Note</b>	For research use only
<b>Formula</b>	$C_{26}H_{30}ClN_7O_2$
<b>MW</b>	508.0151
<b>CAS Number</b>	[1805787-93-2]
<b>SMILES</b>	<chem>O=C(C=C)NC1=CC(CNC2=NC(NC3=C(OC)C=C(N4CCN(C)CC4)C=C3)=NC=C2C)=CC=C1</chem>
<b>Expiration Date</b>	12 months from date of receipt.