

Product Datasheet

Dabigatran ethyl ester (orb611855)

Description	A reversible and selective, direct thrombin inhibitor (DTI) with Ki value of 4.5 nM.
Note	For research use only
Formula	C ₂₇ H ₂₉ N ₇ O ₃
MW	499.5643
CAS Number	[429658-95-7]
SMILES	<chem>C(C1C=CC(NCC2N(C)C3C=CC(C(N(CCC(OCC)=O)C4C=CC=CN=4)=O)=CC=3N=2)=CC=1)(=N)N</chem>
Expiration Date	12 months from date of receipt.