

## Product Datasheet

### (S)-Willardiine (orb611874)

<b>Description</b>	(S)-Willardiine is a potent agonist of AMPA/kainate receptors with EC50 of 44.8 uM..
<b>Note</b>	For research use only
<b>Formula</b>	$C_7H_9N_3O_4$
<b>Purity</b>	> 98%,Standard References
<b>MW</b>	199.1641
<b>CAS Number</b>	[21416-43-3]
<b>SMILES</b>	<chem>NOC(CCN1C=CC(=O)NC1=O)=O</chem>
<b>Expiration Date</b>	12 months from date of receipt.