

Product Datasheet

MLS-1547 (orb611972)

Description	MLS-1547 is a G protein-biased agonist of dopamine D2 receptors (EC50 = 0.37 μ M in a calcium mobili...
Note	For research use only
Formula	$C_{19}H_{19}ClN_4O$
MW	354.8
CAS Number	[315698-36-3]
SMILES	<chem>OC1=C(C=C(Cl)C2=C1N=CC=C2)CN(CC3)CCN3C4=CC=CC=N4</chem>
Expiration Date	12 months from date of receipt.