

## Product Datasheet

### Caspase-1/4 Chromogenic Substrate (orb611985)

<b>Description</b>	Suc-YVAD-pNA is a colorimetric substrate for caspase-1/interleukin-1 $\beta$ -converting enzyme (ICE) and c...
<b>Note</b>	For research use only
<b>Formula</b>	C <sub>31</sub> H <sub>38</sub> N <sub>6</sub> O <sub>12</sub>
<b>MW</b>	686.7
<b>CAS Number</b>	[208264-84-0]
<b>SMILES</b>	<chem>O=C([C@H](CC(O)=O)NC([C@H](C)NC([C@H](C(C)C)NC([C@@H](NC(CC(O)=O)=O)CC1=CC=C(O)C=C1)=O)=O)NC2=CC=C([N+](O-))=O)C=C2</chem>
<b>Expiration Date</b>	12 months from date of receipt.