

Product Datasheet

gly-arg p-nitroanilide dihydrochloride (orb611990)

Description	H-Gly-Arg-pNA is a colorimetric substrate for thrombin.
Note	For research use only
Formula	$C_{14}H_{21}N_7O_4 \cdot 2HCL$
MW	424.28
CAS Number	[125455-61-0]
SMILES	<chem>Cl.NC(NCCCC(C(NC1=CC=C([N+])([O-])=O)C=C1)=O)NC(CN)=O)=N</chem>
Expiration Date	12 months from date of receipt.